



# **Freelite update and Hevylite unpublished data : clinical applications for multiple myeloma**

*Alison Levoguer MB, BS, D.Phil (Oxon)*

*Scientific Affairs Manager, The Binding Site*

# Included in this talk

## Screening for monoclonal gammopathies

Replacement of urine BJP tests at diagnosis (Katzmann, Holding, Cockwell AKI algorithm)


## MGUS prognosis

Risk stratification & 2010 MGUS guidelines (Kyle)


## Monitoring

Multiple myeloma -LCMM (Avet L'Oiseau ASH 2011)

# Monoclonal Gammopathies




# Detection of Monoclonal Proteins



g/L

Monoclonal  
Intact Immunoglobulins

Serum  
electrophoresis



mg/L

Monoclonal  
Free Light Chains

Urine  
electrophoresis  
*or*  
Serum FLC

# Replacement of urine tests with Freelite

Katzmann, J.A., et al., *Clin Chem*, 55: 8; 2009

<b>1877 B cell disorders</b>	
Multiple myeloma	467
AL amyloidosis	581
Light chain deposition disease	18
Waldenström's macroglobulinemia	26
Plasmacytoma	29
Extramedullary plasmacytoma	10
POEMS	31
Smouldering myeloma	191
MGUS	524

# Replacement of urine tests with Freelite

Screening algorithm	Diagnostic sensitivity (%)		Negative screens compared to all 5 tests
SPE, sIFE, UPE, uIFE sFLC	MM AL LCDD WM	100 98.1 83.3 100	

# Replacement of urine tests with Freelite

Screening algorithm	Diagnostic sensitivity (%)		Negative screens compared to all 5 tests
<b>SPE, sIFE, UPE, uIFE sFLC</b>	<b>MM</b>	100	
	<b>AL</b>	98.1	
	<b>LCDD</b>	83.3	
	<b>WM</b>	100	
<b>SPE alone</b>	<b>MM</b>	87.6	58
<b>+ reflex sIFE</b>	<b>AL</b>	67.2	187
	<b>LCDD</b>	66.7	5
	<b>WM</b>	100	0

# Replacement of urine tests with Freelite

Screening algorithm	Diagnostic sensitivity (%)		Negative screens compared to all 5 tests
<b>SPE, sIFE, UPE, uIFE sFLC</b>	<b>MM</b>	100	
	<b>AL</b>	98.1	
	<b>LCDD</b>	83.3	
	<b>WM</b>	100	
<b>SPE alone + reflex sIFE</b>	<b>MM</b>	87.6	58
	<b>AL</b>	67.2	187
	<b>LCDD</b>	66.7	5
	<b>WM</b>	100	0
<b>SPE, sIFE + uIFE</b>	<b>MM</b>	98.7	6
	<b>AL</b>	94.2	23
	<b>LCDD</b>	77.8	1
	<b>WM</b>	100	0



# Replacement of urine tests with Freelite

Screening algorithm	Diagnostic sensitivity (%)		Negative screens compared to all 5 tests
<b>SPE, sIFE, UPE, uIFE sFLC</b>	<b>MM</b>	100	
	<b>AL</b>	98.1	
	<b>LCDD</b>	83.3	
	<b>WM</b>	100	
<b>SPE alone + reflex sIFE</b>	<b>MM</b>	87.6	58
	<b>AL</b>	67.2	187
	<b>LCDD</b>	66.7	5
	<b>WM</b>	100	0
<b>SPE, sIFE + uIFE</b>	<b>MM</b>	98.7	6
	<b>AL</b>	94.2	23
	<b>LCDD</b>	77.8	1
	<b>WM</b>	100	0
<b>SPE + sFLC + reflex sIFE</b>	<b>MM</b>	100	0
	<b>AL</b>	96.2	11
	<b>LCDD</b>	77.8	1
	<b>WM</b>	100	0

# Does screening with Freelite have any additional value?

In the screening algorithm of a general hospital

## SPE/CZE+Urine vs SPE/CZE+Freelite

**Hill *et al. Clin Chem* (2006); 52: 1743-1748**

Detection increased by 12%

**Holding *et al. Clin Chem Lab Med* (2011); 49:83-88**

Detection increased by 7%

## SPE/CZE vs SPE/CZE+Freelite

**Bakshi *et al. Am J Clin Pathol* (2005); 124: 214-218**

Detection increased by 41%

**Robson *et al. Lab Med* (2009); 40: 325-329**

Detection increased by 71%

# UK Hull screening study

Prospective study:

753 samples received for MG screen

FLC assays detected 5 additional pts negative by SPE/CZE (2 MM and 3 MGUS)

For these 5 pts, only 3 urine samples received (2, 30 and 49 days after serum)

Overall, only 17% of SPE requests had UPE performed within 90 days

# Acute Kidney Injury (AKI)

Screening algorithm

**Unexplained  
acute kidney injury**

**Screen by sFLC**

**Clonal FLC  
> 500 mg/L**

**Probable  
CN**

**Clonal FLC  
< 500 mg/L**

**No  
CN**

**Normal  
FLC ratio**

**Monoclonal  
pathology  
unlikely**

**Consider renal biopsy  
Haemato-oncology work-up**

# Included in this talk

## Screening for monoclonal gammopathies

Replacement of urine BJP tests at diagnosis (Katzmann, Holding, Cockwell AKI algorithm)

## MGUS prognosis

Risk stratification & 2010 MGUS guidelines (Kyle)

## Monitoring

Multiple myeloma -LCMM (Avet L'Oiseau ASH 2011)

# Serum FLCs are prognostic in:

- MGUS progression
- Smouldering MM progression
- Plasmacytoma progression
- Myeloma outcome
- AL amyloidosis outcome
- B-CLL outcome
- Waldenström's outcome

# MGUS risk stratification model incorporating M-protein size, type and FLC ratio

Risk of progression	No. of abnormal risk factors	No. patients	Absolute risk of progression at 20 years*
Low	0	449	2%
Low-Intermediate	1	420	10%
High-Intermediate	2	226	18%
High	3	53	27%

\* Accounting for death as a competing risk

**Rajkumar et al. *Blood* 2005: 106, 812-817**




# 2010 IMWG MGUS Guidelines

Risk of progression	Recommended follow-up
Low	6 months then every 2-3 years <u>or at progression</u>
Low-Intermediate	6 months then annually
High-Intermediate	
High	


**Kyle et al. *Leukemia* 2010: 24, 1121-1127**

# Monitoring Patients Following Treatment


*Which patients can be monitored with Freelite following treatment?*




Nonsecretory MM



Light chain only MM



AL Amyloidosis




Intact Immunoglobulin MM

89% of IIMM also have abnormal FLC ratio

Mead et al. Br J Haematol 2004;126 : 348 – 54

# Monitoring LCMM



# Summary

- sFLC assays can replace urine testing in the diagnosis of monoclonal gammopathies
- sFLC assays improve screening of AKI
- sFLC ratio is an independent risk factor for MGUS progression to malignant disease
- sFLC assays assist with management as they
  - rapidly indicate response to treatment
  - allow detection of light chain escape




## Hevylite recent data

# Focus of this part of the talk




- An overview of HLC analysis & its benefits
- Published & unpublished clinical data
  - Monitoring
  - Prognosis

# Detection of Monoclonal Proteins



Monoclonal  
Intact Immunoglobulins

Serum electrophoresis  
or  
Serum HLC



Monoclonal  
Free Light Chains

Urine electrophoresis  
or  
Serum FLC

# IMWG guidelines: Multiple myeloma follow-up

- “For patients with measurable monoclonal protein in serum, both electrophoretic studies and quantitative immunoglobulins are recommended to assess response, although electrophoretic measurements to follow monoclonal protein are preferred.”




# Problems with electrophoresis

- Some serum protein electrophoresis (SPE) results are difficult to interpret
  - Broad bands & smearing of bands
  - Couplet/ triplets of bands
  - Precipitation at the point of application
  - Dye saturation
  - Hidden bands

# Problems with SPE densitometry

- IgA bands may be hidden by transferrin



# IMWG guidelines:

## Multiple myeloma follow-up

- “For patients with measurable monoclonal protein in serum, both electrophoretic studies and quantitative immunoglobulins are recommended to assess response, although electrophoretic measurements to follow monoclonal protein are preferred.”
- “For several patients, especially with IgA or IgD myeloma, nephelometric quantitation of serum immunoglobulin is necessary.”

# Standard nephelometric IgG, IgA and IgM assays


Total immunoglobulin measurements do not distinguish monoclonal from polyclonal proteins within any given sample

e.g. Total IgA =

$$\text{IgA}\kappa + \text{IgA}\lambda$$

*involved* immunoglobulin + *uninvolved* immunoglobulin

# Introducing Hevylite analysis




- Polyclonal sheep antisera reagents
- Specific for junctional epitopes that span immunoglobulin heavy & light chains

IgA $\kappa$

IgA $\lambda$

# Hevylite specificities




IgG $\kappa$     IgG $\lambda$



IgG $\kappa$ / IgG $\lambda$

NR=1.12-3.21\*




IgA $\kappa$     IgA $\lambda$



IgA $\kappa$ / IgA $\lambda$

NR= 0.78-1.94\*



IgM $\kappa$     IgM $\lambda$



IgM $\kappa$ / IgM $\lambda$

NR=1.18-2.74\*

\* Generated on SPAPLUS analyser

# Additional information from HLC measurements

Involved immunoglobulin & uninvolved immunoglobulin

A large, bright green circle containing the text "IgAκ" in a large, black, serif font. The text is centered within the circle.

A smaller, medium-green circle containing the text "IgAλ" in a black, serif font. Below the text, the phrase "HLC pair suppression" is written in a smaller, black, sans-serif font. The text is centered within the circle.


# Another potential advantage of IgG HLC ratios

## Corrects for variable IgG metabolism


- IgA 6 days
- IgM 5 days
- **IgG** ~ **20 days** (*days to weeks*)



# FcRn binds IgG




# Concentration-dependent catabolism of IgG




*HLC ratio corrects for this*

# IgA HLC at presentation



# IgG HLC at presentation



Using HLC assays in routine practice

Monitoring multiple myeloma

# A 71-Year-Old Woman with Multiple Myeloma Status after Stem Cell Transplantation


Leslie J. Donato,<sup>1</sup> Steven R. Zeldenrust,<sup>2</sup> David L. Murray,<sup>1</sup> and Jerry A. Katzmann<sup>1,2\*</sup>

## Clinical History

- 9 yr history of MGUS

## Presentation

- Anaemia and hypercalcaemia
- SPE - 48 g/L M protein
- IFE – monoclonal IgAk
- Raised total IgA (47.2g/L)
- Abnormal sFLC  $\kappa/\lambda$  ratio (7)
- BM - 40% plasma cells



# A 71-Year-Old Woman with Multiple Myeloma Status after Stem Cell Transplantation

Leslie J. Donato,<sup>1</sup> Steven R. Zeldenrust,<sup>2</sup> David L. Murray,<sup>1</sup> and Jerry A. Katzmann<sup>1,2\*</sup>

## Diagnosis


- MM (Durie-Salmon stg IIIA/Int stg 2)

## Treatment & monitoring

- Stem cell transplant
- Pt became asymptomatic
- Electrophoresis negative for 1.5 yrs

## 1<sup>st</sup> relapse

- 2 yrs following ASCT




**Laboratory results for IFE interpretation, immunoglobulin FLC ratio, IgA  
HLC ratio and IgA patient monitoring**  
Clinical Chemistry 57:12 (2011) 1645-9


Timepoint	IgA g/L	IFE	FLC ratio	Hevylite IgAk/IgAl ratio	Disease status
Reference interval	.610-3.56	-ve	0.26 – 1.65	0.7 – 2.2	
1	38.1	M-spike 39.2 g/L	10.3	463	1 <sup>st</sup> relapse
2	3.28	-ve	1.9	3.7	Rx change
3	4.34	-ve	1.14	2.9	Side effects
4	4.28	-ve	1.06	2.8	Revlimid
5	6.58	-ve	1.21	4.1	Revlimid
6	7.67	-ve	1.08	5.6	Revlimid
7	9.14	-ve	1.36	5.6	+ Dex
8	8.04	-ve	1.38	8.2	Progression
9	8.80	+ve	2.07	13.0	Relapse
10	12.20	+ve	4.75	24.8	Rx change



# Monitoring over 3 years (Jan 2008 – Jan 2011)



# IgAκ/IgAλ HLC ratio NEVER normalised



---

# A 71-Year-Old Woman with Multiple Myeloma Status after Stem Cell Transplantation

Leslie J. Donato,<sup>1</sup> Steven R. Zeldenrust,<sup>2</sup> David L. Murray,<sup>1</sup> and Jerry A. Katzmann<sup>1,2\*</sup>

---

“Monoclonal proteins that migrate as broad bands can be difficult to distinguish, and quantitative assessment of immunoglobulin HLC pairs provides a measure of clonal synthesis.”

Donato et al. Clin Chem. 2011;57(12):1645-8

**Commentary** - James D. Faix, Stanford University

“Now it is clear that we have another alternative for cases like this.”

**PRACA ORYGINALNA – Original Article**

MARIA KRAJ, BARBARA KRUK, RYSZARD POGLÓD, KRZYSZTOF WARZOCHA

**Evaluation of IgG, IgA and IgM monoclonal and biclonal gammopathies by nephelometric measurement of individual immunoglobulin  $\kappa/\lambda$  ratios – Hevylite assay versus immunofixation**

**Ocena IgG, IgA i IgM monoklonalnej i biklonalnej gammopatii metodą nefelometrycznego oznaczania stosunku  $\kappa/\lambda$  indywidualnych immunoglobulin – Test Hevylite versus immunofiksacja**


Institute of Hematology and Transfusion Medicine, Warsaw, Poland  
Head: Prof. Krzysztof Warzocha

---

Why FLC + HLC in MM patients?


# Why FLC and HLC?

sFLC and intact immunoglobulins are independent tumour markers.....



...and may be produced by different plasma cell clones


IgG $\kappa$  multiple myeloma:



IgG $\kappa$


$\kappa$  FLC

**Freelite** and **Hevylite** to monitor all  
plasma cell tumour clones






# Why FLC and HLC?




# Using HLC assays in routine practice

Early prognostic data

# Negative correlation between involved and uninvolved IgG HLC




# Extreme Hevylite ratios



# Prognostic value of HLC analysis

Baseline HLC ratios (IgG  
and IgA MM)




# 3-tier risk stratification model

## Risk Factors

HLC ratio  $<0.022$  or  $>45$

FLC ratio  $<0.1$  or  $>30$



n = 310

5 year survival rates:

**0 risk factors = 67.4%**

**1 risk factor = 50.0%**

**2 risk factors = 23.3%**

# Hevylite: key points



- Simple, fast, automated tests
- Clear identification, quantitation and typing of M - protein
- HLC ratio is a sensitive measure of disease and takes into account both the monoclonal and underlying polyclonal levels
- Early data indicates that the HLC ratio is prognostic

# Are you interested in a clinical study?

- Hard-to-monitor patients?
- Monitoring protocol can be provided
- Prognostic benefit





- Chapter 25 – Guidelines for use of sFLC assays
  - 25.5. IMWG: guidelines for standard investigative workup of patients with suspected multiple myeloma (2011)
  - 25.6. IMWG: guidelines for risk stratification in multiple myeloma (2011)
  - 25.7. IMWG: consensus recommendations for the uniform reporting of clinical trials (2011)
  
- Chapter 32 – Analysis of Hevylite

# Any questions?



- [alison.levoguer@bindingsite.co.uk](mailto:alison.levoguer@bindingsite.co.uk)

**Keep up-to-date at:**

**Wikilite**  
*.COM*





# IgA HLC: Monitoring IgA MM




*IgA HLC ratio has **greater sensitivity** than IFE for detection of residual disease in some IgA patients*

# IgG HLC: Monitoring IgG MM



- IgG Hevylite ratio is normal when IFE is positive for some IgG patients
  - ?Hevylite insensitive due to background IgG
  - ?IFE is false positive due to recycling by FcRN receptor

# IFE positivity remains due to IgG recycling



*More studies are needed to determine the clinical significance*

# Hevylite & biclonal gammopathies

- Bicloneals are rare (approx 1%-5% of clonal gammopathies)
- For biclone of different isotypes (ie IgGK & IgAK) use both Hevylite pairs for monitoring
- Hevylite ratio may be normal for biclone of same isotype (ie IgGK & IgG $\lambda$ ) so use SPE/IFE for monitoring these rare cases


# • Hevylite + Freelite are

- Independent biomarkers
- Used together both give different clinical information
- Run both together to optimize clinical practice




# Freelite & Hevylite: Patient monitoring


	Freelite	Hevylite
Nonsecretory/ Oligosecretory MM		?
Light chain MM		?
Intact Immunoglobulin MM		




Intact IgGk



κ FLC



γ heavy chain









Binding site for Hevylite reagent

Binding site for Freelite reagent

κ  
Freelite  
binding?

IgGκ  
Hevylite  
binding?

 Epitopes hidden by bound heavy chain	 Epitopes present
 Epitopes now exposed	 Epitopes absent
 Epitopes absent	 Epitopes absent

# Limitations of total IgG, IgA and IgM assays

Consider a total IgG result of 15g/L

9g/L IgG $\kappa$  + 6g/L IgG $\lambda$

14.5g/L IgG $\kappa$  + 0.5 g/L IgG $\lambda$


IgG $\kappa$ / IgG $\lambda$

1.5

29

A ratio rather than total value is much more informative...


# Monitoring with IgA Hevylite



# Normal Ratio Ranges

Assays	95% range (g/L)
IgG $\kappa$ / IgG $\lambda$ ratio	1.12-3.21
IgA $\kappa$ / IgA $\lambda$ ratio	0.78-1.94
IgM $\kappa$ / IgM $\lambda$ ratio	1.18-2.74

# Recycling of IgG



Courtesy of J Hobbs